

Water Testing



AquaScreen®

DNA-based system for quantitative detection of water pathogens.

AquaScreen® combines water filtration, lysis of the collected microorganisms, DNA extraction and elution of the DNA in minimal volumes ready for PCR analysis.

Very sensitive

Detection limit comparable to culture method. Sample volume up to 1000 ml applicable.

Ease of use

Sample preparation includes simple water filtration and DNA extraction prior to direct PCR analysis.

Instrument Compatibility

Our detection kits can be used on most commercially available qPCR cyclers equipped with standard filters for FAM™ and ROX™.

Label Statement

Fully licensed probe system.

One extract – multiple parameters

The AquaScreen® FastExtract procedure obtains suitable template DNA for routine PCR tests. One extract can be used for PCR reactions with different specificities, so that multiple microorganisms can be analyzed in parallel and on users' choice.

Fast

The fast extraction procedure in combination with the rapid PCR method allows detection of pathogenic organisms (including those that may be non-culturable) in industrial and environmental water samples within a few hours.

Recognizes VBNC (Viable But Not Culturable)

Free DNA of dead and already lysed microorganisms in the sample passes through the membrane filter and is not detectable in the test system. The test detects both "viable and culturable" and "viable but not culturable" (VBNC) legionella. The effect of hygiene measures like heat treatment can be monitored after appropriate flushing and recovery time.

Clear

A clear and easily interpretable result is obtained with one PCR reaction. No subsequent and laborious specification methods are required.

Stable

All PCR kits are Freeze-dried and need to be rehydrated with a supplied buffer to reduce shipping costs and increase product stability.

Products for the control of microbials in
cell cultures & water.



Water Testing

Procedure

Step 1: Preparation of the Sample Material

Drinking water, condensed water from cooling systems, bathing or pool water, and waste water released from suspended particles can be used as sample material. As intact microorganisms are needed for the filtration procedure thermal decontamination procedures can not be monitored immediately. With suspended particles or fixed volatile contents, contaminated water samples can be purified by prior filtration with a paper folded filter. The samples may not be centrifuged for purification. For the testing procedure, at least 100 ml is minimally required, however a sample volume of 1000 ml is recommended for highest sensitivity.

Step 2: PCR Application

Currently AquaScreen® qPCR kits are available for quantitative detection of *Legionella pneumophila*, *Legionella species*, *Pseudomonas aeruginosa* and *Escherichia coli* in water samples. Kits for additional parameter are in preparation. The test is based on quantitative real-time PCR (qPCR) which allows for highest specificity and sensitivity. The PCR mix contains a primer/probe set specific for the microorganism to be detected and emits a fluorescent light at ~520 nm. The kit includes an internal control that is detected by another probe emitting light at ~560 nm. The internal control is constructed as a homologous control containing primer binding sites identical to the target but with a unique internal sequence. By using the supplied internal control, false negative results (e.g. due to inhibition of the reaction by the sample matrix) can be excluded individually for each sample.

AquaScreen® FastExtract



For latest information on our website please scan the QR code



Features		Features	
Description	Rapid DNA extraction from water samples	Package Sizes	Cat.-No. 32-1010 10 extractions Cat.-No. 32-1050 50 extractions Cat.-No. 32-1200 4x50 extractions
Recommended Use / Scope	AquaScreen® FastExtract can be used with your established suction device (47 mm frit) for the extraction of legionella and other microbial contaminations. AquaScreen® FastExtract is optimized for high flow and throughput and provides high quality DNA for subsequent PCR analysis.	Required lab devices & reagents	Vacuum pump Micro centrifuge Filtration system, 47 mm frit Pipetting equipment and filtered tips Incubator (37 °C for petri dishes, 56 °C for reaction tubes) Ethanol (96-100 %)
Kit Components	Membrane filters Incubation dishes Incubation, collection and sample storage tubes Lysis, wash and elution buffers	Shelf Life and Storage	Components are maintainable at room temperature for at least 6 months.
		Compliance	AFNOR XP T90-471 and ISO/TS 12869:2012 in combination with AquaScreen® qPCR kits

AquaScreen® Legionella pneumophila



For latest information on our website please scan the QR code



Features		Features	
Type of PCR	Quantitative real-time PCR (qPCR)	Package Sizes	Cat.-No. 34-2025 25 reactions Cat.-No. 34-2100 100 reactions Cat.-No. 34-2250 250 reactions
Description	The AquaScreen® Legionella pneumophila qPCR Detection Kit is used for DNA samples prepared with the AquaScreen® FastExtract procedure for quantification of <i>Legionella pneumophila</i> in waters samples. The supplied primer set is specific for a segment of the mip region of the <i>Legionella pneumophila</i> genome.	Required Consumables	PCR reaction tubes Optional: For calibration we recommend our <i>Legionella pneumophila</i> DNA Calibration Reagent (Cat.-No. 52-0101).
Recommended Use / Scope	Applicable in research and industry for QA testing of household and process water. Not recommended for clinical diagnostics, testing of human samples or pharmaceutical products.	Required lab devices	Pipetting equipment qPCR cycler with filter sets for FAM™ and ROX™
Kit Components	Freeze-dried primer, probes, nucleotides and polymerase Rehydration Buffer Freeze-dried Positive Control DNA Freeze-dried Internal Amplification Control	Shelf Life and Storage	Components are maintainable at +2 to +8 °C for at least 6 months. After rehydratisation the reagents must be stored at -18 °C
		Compliance	AFNOR XP T90-471 and ISO/TS 12869:2012

AquaScreen® Legionella species



For latest information on our website please scan the QR code



Features		Features	
Type of PCR	Quantitative real-time PCR (qPCR)	Package Sizes	Cat.-No. 33-2025 25 reactions Cat.-No. 33-2100 100 reactions Cat.-No. 33-2250 250 reactions
Description	The AquaScreen® Legionella species qPCR Detection Kit is used for DNA samples prepared with the AquaScreen® Fast Extract procedure for quantification of <i>Legionella</i> in waters samples. The supplied primer set is specific for a broad range of legionella species, but does not detect other water born bacteria as required by ISO/TS 12869:2012.	Required Consumables	PCR reaction tubes Optional: For calibration we recommend our <i>Legionella pneumophila</i> DNA Calibration Reagent (Cat.-No. 52-0101).
Recommended Use / Scope	Applicable for water testing as described in ISO/TS 12869:2012, in research and industry for QA testing of process water. Not recommended for clinical diagnostics, testing of human samples or pharmaceutical products.	Required lab devices	Pipetting equipment qPCR cycler with filter sets for FAM™ and ROX™
Kit Components	Freeze-dried primer, probes, nucleotides and polymerase Rehydration Buffer Freeze-dried Positive Control DNA Freeze-dried Internal Amplification Control	Shelf Life and Storage	Components are maintainable at +2 to +8 °C for at least 6 months. After rehydratisation the reagents must be stored at -18 °C
		Compliance	AFNOR XP T90-471 and ISO/TS 12869:2012

Water Testing

AquaScreen® Pseudomonas aeruginosa



For latest information on our website please scan the QR code



Features		Features	
Type of PCR	Quantitative real-time PCR (qPCR)	Package Sizes	Cat.-No. 34-6025 25 reactions Cat.-No. 34-6100 100 reactions Cat.-No. 34-6250 250 reactions Primer sets and nucleotides are prepared in aliquots of 25 tests.
Description	The AquaScreen® Pseudomonas aeruginosa qPCR Detection Kit is used for DNA samples prepared with the AquaScreen® FastExtract procedure for quantification of Pseudomonas aeruginosa in waters samples.	Required Consumables	PCR reaction tubes Optional: For calibration we recommend our Pseudomonas aeruginosa DNA Calibration Reagent (Cat.-No. 52-0071).
Recommended Use / Scope	Applicable in research and industry for QA testing of household and process water. Not recommended for clinical diagnostics, testing of human samples or pharmaceutical products.	Required lab devices	Pipetting equipment qPCR cyclers with filter sets for FAM™ and ROX™
Kit Components	Freeze-dried primer, probes, nucleotides and polymerase Rehydration Buffer Freeze-dried Positive Control DNA Freeze-dried Internal Amplification Control	Shelf Life and Storage	Components are maintainable at +2 to +8 °C for at least 6 months. After rehydration the reagents must be stored at -18 °C.
		Compliance	No guidelines are available for molecular testing of water samples for Pseudomonas aeruginosa.

AquaScreen® Escherichia coli



For latest information on our website please scan the QR code



Features		Features	
Type of PCR	Quantitative real-time PCR (qPCR)	Package Sizes	Cat.-No. 34-7025 25 reactions Cat.-No. 34-7100 100 reactions Cat.-No. 34-7250 250 reactions Primer sets and nucleotides are prepared in aliquots of 25 tests.
Description	The AquaScreen® Escherichia coli qPCR Detection Kit is used for DNA samples prepared with the AquaScreen® FastExtract procedure for quantification of Escherichia coli in waters samples.	Required Consumables	PCR reaction tubes Optional: For calibration we recommend our Escherichia coli DNA Calibration Reagent (Cat.-No. 52-0083).
Recommended Use / Scope	Applicable in research and industry for QA testing of household and process water. Not recommended for clinical diagnostics, testing of human samples or pharmaceutical products.	Required lab devices	Pipetting equipment qPCR cyclers with filter sets for FAM™ and ROX™
Kit Components	Freeze-dried primer, probes, nucleotides and polymerase Rehydration Buffer Freeze-dried Positive Control DNA Freeze-dried Internal Amplification Control	Shelf Life and Storage	Components are maintainable at +2 to +8 °C for at least 6 months. After rehydration the reagents must be stored at -18 °C.
		Compliance	No guidelines are available for molecular testing of water samples for Escherichia coli.