

# 5

## reasons

that make these tubes...

## Simply Perfect!

- **Easy seal design** for repeated single hand opening & closing. Smooth cap won't catch gloves.
- **Clear large graduations** from 0.1ml to 1.5ml on the body.
- Very clear **transparent surface**. Doesn't hinder spectrophotometer reading.
- Bottom designed with **no sharp points** to avoid glove tearing.
- **Thin cap membrane** for trouble-free needle piercing.



Manufactured from FDA approved resins  
Certified RNase, DNase, DNA & Pyrogen free  
Can be centrifuged at upto 20,000 RCF  
Withstands temperatures from -80°C to +121°C  
Fit all standard and high capacity centrifuge rotors  
Supplied in double safety self standing bulk bags



ISO 9001:2008  
ISO 13485:2003

Cat. No.	Description	Packing	Qty / Pack	Qty / Case
ACT-017-B	1.7 ml Graduated Micro Centrifuge Tube	Bulk	1000 Tubes	1000 Tubes X 10 packs



Accumax Lab Technology

[www.accumaximum.com](http://www.accumaximum.com)