

# Accumax

## Pipette Tips

Produced from high quality diamond finish mould using special medical grade virgin polypropylene with very high transparency

Innovative design of tips for adaptability with wide range of pipettes with minimum force requirement for loading & ejection

Ideal match for high sensitive PCR work

Produced in a fully automatic, human touch free, clean room production facility. Free from DNase, RNase, Human DNA, Pyrogen & PCR inhibitor. Free of Lubricants and heavy metals

Graduation marks on all the tips that act as visual checks to minimize any risk of errors during sample aspiration

All these high quality tips are also available as sterile tips



10 µl Universal Graduated Tips				
1000 Tips / Bag (40 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-10-B	-	-	-	
AT-10-S-B	✓	-	-	
AT-10-L-B	-	✓	-	
AT-10-F-B	-	-	✓	
AT-10-S-L-B	✓	✓	-	
AT-10-S-F-B	✓	-	✓	
AT-10-L-F-B	-	✓	✓	
AT-10-S-L-F-B	✓	✓	✓	

20 µl Universal Graduated Tips				
(With Filter only)				
1000 Tips / Bag (20 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-20-F-B	-	-	✓	
AT-20-S-F-B	✓	-	✓	
AT-20-L-F-B	-	✓	✓	
AT-20-S-L-F-B	✓	✓	✓	

50 µl Universal Graduated Tips				
(With Filter only)				
1000 Tips / Bag (20 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-50-F-B	-	-	✓	
AT-50-S-F-B	✓	-	✓	
AT-50-L-F-B	-	✓	✓	
AT-50-S-L-F-B	✓	✓	✓	

## BULK


Packed in polyethylene bag with easy tearing

**Perforation & Zip lock**- Double sealing provides better protection against contamination and spillage while in use

100 µl Universal Graduated Tips				
(With Filter only)				
1000 Tips / Bag (20 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-100-F-B	-	-	✓	
AT-100-S-F-B	✓	-	✓	
AT-100-L-F-B	-	✓	✓	
AT-100-S-L-F-B	✓	✓	✓	

200 µl Universal Graduated Tips				
1000 Tips / Bag (20 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-200-B	-	-	-	
AT-200-S-B	✓	-	-	
AT-200-L-B	-	✓	-	
AT-200-F-B	-	-	✓	
AT-200-S-L-B	✓	✓	-	
AT-200-S-F-B	✓	-	✓	
AT-200-L-F-B	-	✓	✓	
AT-200-S-L-F-B	✓	✓	✓	

1000 µl Universal Graduated Tips				
1000 Tips / Bag (10 Bags / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-1000-B	-	-	-	
AT-1000-S-B	✓	-	-	
AT-1000-L-B	-	✓	-	
AT-1000-F-B	-	-	✓	
AT-1000-S-L-B	✓	✓	-	
AT-1000-S-F-B	✓	-	✓	
AT-1000-L-F-B	-	✓	✓	
AT-1000-S-L-F-B	✓	✓	✓	



## RACK

Fully autoclavable and stackable

Hinge rack allows single hand operation

Color coded tray for easy identification

10 µl Universal Graduated Tips				
96 Tips x 10 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-10-R	-	-	-	
AT-10-S-R	✓	-	-	
AT-10-L-R	-	✓	-	
AT-10-F-R	-	-	✓	
AT-10-S-L-R	✓	✓	-	
AT-10-S-F-R	✓	-	✓	
AT-10-L-F-R	-	✓	✓	
AT-10-S-L-F-R	✓	✓	✓	

20 µl Universal Graduated Tips				
(With Filter only)				
96 Tips x 10 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-20-F-R	-	-	✓	
AT-20-S-F-R	✓	-	✓	
AT-20-L-F-R	-	✓	✓	
AT-20-S-L-F-R	✓	✓	✓	

50 µl Universal Graduated Tips				
(With Filter only)				
96 Tips x 10 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-50-F-R	-	-	✓	
AT-50-S-F-R	✓	-	✓	
AT-50-L-F-R	-	✓	✓	
AT-50-S-L-F-R	✓	✓	✓	

100 µl Universal Graduated Tips				
(With Filter only)				
96 Tips x 10 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-100-F-R	-	-	✓	
AT-100-S-F-R	✓	-	✓	
AT-100-L-F-R	-	✓	✓	
AT-100-S-L-F-R	✓	✓	✓	

200 µl Universal Graduated Tips				
96 Tips x 10 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-200-R	-	-	-	
AT-200-S-R	✓	-	-	
AT-200-L-R	-	✓	-	
AT-200-F-R	-	-	✓	
AT-200-S-L-R	✓	✓	-	
AT-200-S-F-R	✓	-	✓	
AT-200-L-F-R	-	✓	✓	
AT-200-S-L-F-R	✓	✓	✓	

1000 µl Universal Graduated Tips				
96 Tips x 6 Racks / Pack (10 Packs / Case)				
Cat #	Sterile	Low Retention	Filter	
AT-1000-R	-	-	-	
AT-1000-S-R	✓	-	-	
AT-1000-L-R	-	✓	-	
AT-1000-F-R	-	-	✓	
AT-1000-S-L-R	✓	✓	-	
AT-1000-S-F-R	✓	-	✓	
AT-1000-L-F-R	-	✓	✓	
AT-1000-S-L-F-R	✓	✓	✓	

DNase

Rnase

Human DNA

Endotoxin

## FILTER TIPS



Specially developed PE filters stops aerosols & contamination from entering into the pipette and thus lessens the risk or cross contamination.

Free from chemical additives & heavy metals

## LOW RETENTION TIPS



Regular tips

Accumax Low Retention tips

**Minimum retention = Maximum accuracy**  
**No silicon additives = No leaching**

Manufactured from specifically developed resins combined with high quality diamond finish mould creates **SUPER HYDROPHOBIC** internal surface

## RELOADING STACK

Reloading of rack with minimal effort

Easy to operate without the need of any special equipment

Tips prearranged in 96 (8x12) format

Saves storage space and minimize packaging waste

### 10 µl Universal Graduated Tips

96 Tips x 10 Plates / Pack (10 Packs / Case)

Cat #	Sterile	Low Retention
AT-10-K	-	-
AT-10-S-K	✓	-
AT-10-L-K	-	✓
AT-10-S-L-K	✓	✓

### 200 µl Universal Graduated Tips


96 Tips x 10 Plates / Pack (10 Packs / Case)

Cat #	Sterile	Low Retention
AT-200-K	-	-
AT-200-S-K	✓	-
AT-200-L-K	-	✓
AT-200-S-L-K	✓	✓

### 1000 µl Universal Graduated Tips


96 Tips x 5 Plates / Pack (10 Packs / Case)

Cat #	Sterile	Low Retention
AT-1000-K	-	-
AT-1000-S-K	✓	-
AT-1000-L-K	-	✓
AT-1000-S-L-K	✓	✓



#### Special Note

Every version of 200µl and 1000µl tips is available in Yellow and Blue color respectively on request. Suffix "C" in catalog# while placing order



ISO 9001:2008

ISO 13485:2003

Accumax Lab Technologies

info@accumaximum.com  
 www.accumaximum.com

*Accumax*

Pipette Tips